



ICE FLAKING MACHINE



SALIENT FEATURES

- Fully automatic microprocessor control
 - Compact Design and low maintenance
 - SS construction
 - Continuous Ice Flakes Output
 - Noise level ≤60 dB which results in silent operation
 - CFC FREE ECO FRIENDLY COMPRESSOR.
 - Well balanced refrigeration system with hermetically sealed compressor, condenser, fan motor assembly
 - Low water level detection
 - Water supply using 3/4" BSP
- Ice flakes coming out from the crushing assembly only within 25 minutes and are collected in a collecting bin for ready use
 - Requires plain tap water connection
 - PUF insulation for minimal thermal loss
 - Drain arrangement at the bottom
 - Crusher assembly is made of Stainless Steel and is driven by heavy duty induction motor coupled to reduction gear.
 - Wheels for ease of movement
 - It works on 230V AC.

SPECIFICATION

Model	SAII-IFM-15	SAII-IFM-30	SAII-IFM -50	SAII-IFM-100	SAII-IFM-150	SAII-IFM-300
Capacity (per 24 hrs)	15	30	50	100	150	300
Storage bin	ABS plastic					
Storage capacity (Kg)	4	10	15	25	30	40
Cooling method	Air					
Working temp.	Room temp below 30 °C					
Input power (W)	140	180	280	420	600	1000
Dimension (cm)	27 x 45 x 50	30 x 50 x 55	38 x 55 x 72	55 x 63 x 90	70 x 70 x 110	90 x 90 x 150
Gross weight (Kg)	25	32	45	65	90	160

Optional:

- Compatible Voltage Stabilizer